

CAROSPEED®

Hand-operated
accelerator
and brake
control



 **auto**
ADAPT
for your independence


INVENTIONS FOR QUALITY OF LIFE

CAROSPEED®

- gives you the right feeling

Right driving feeling

Carospeed is a hand-operated accelerator and brake control for drivers who are unable to use the car pedals. The accelerator and brake functions are transmitted to a lever that is placed beside the driver's seat. Pull the lever backwards and the vehicle accelerates, brake by pushing the lever forward. It is so simple to handle. Carospeed gives the driver the right feeling when accelerating and braking. The secret lies in the lever's geometry and our long experience in vehicle adaptations.



Carospeed is a universal hand control that can be fitted in most vehicles, both right and left hand drive, large and small. Carospeed is fixed to the floor using the drivers seat fixtures and is connected to the ordinary accelerator and brake pedals.

Hill-holder as standard

The hand control is equipped with a comfortable hill-holder as standard. The hill-holder keeps the vehicle still while the hand moves between gear lever, regulator or hand brake.

Set function for cruise

To avoid keeping the hand on the accelerator when driving far distances, the hand control is prepared for fitting in conjunction with the vehicle's original cruise control. The set function for the cruise control is activated during driving with the same button as the hill-holder and is disconnected automatically when braking. (The fitting of the set function is an option).

Ergonomic design

Carospeed is an ergonomic designed aid. Many years of experience in designing and manufacturing vehicle adaptations lies behind the development of Carospeed. Demands for easy installation and optimal drive feeling has ruled the development. A reference group of disabled drivers has influenced the development and we have had high demands on the hand control design to match vehicle interior equipment.

CAROSPEED®

Hill-holder and set function for cruise control.



CAROSPEED®



Black handle and grey cover around the lever bar.



Wood/black handle and black cover around the lever bar.

...eling when driving!



Safety first

Carospeed is CE-marked and meets safety requirements from all major car manufacturers. It is also e-marked within the European community (95/54/EC).



Carospeed is accredited in the UK by Motability.



The installation must be carried out by trained technicians from an Autoadapt authorized workshop.



The interference is minimal and the vehicle can easily be restored when the time comes to change. Carospeed is recommended for vehicles with automatic transmission.

CAROSPEED® T

CAROSPEED® TE



Carospeed T (Standard) and **TE** (Electrical functions) are hand controls with a T-grip. The handle is slightly sloping, which gives the hand a restful and ergonomic position.

Carospeed T and TE are always delivered with a choice of two different handle designs. One design is a soft cylinder shape and the other is a modelled design with solid surface. This way the driver can choose which one he/she feels is the most comfortable.

Wind screen washer, interval, wiper 1, wiper 2



Direction indicator

Headlight

Horn

Hill-holder and set function for cruise control.



CAROSPEED® E

CAROSPEED® INDICATOR

Carospeed Indicator has a switch only for the turn left and the turn right indicator.



...dle and grey lever bar.



T-grip with a soft cylinder shaped design and T-grip in a modelled design with solid surface.

Hill holder and set function for cruise are standard. Carospeed Indicator is available with a black coloured handle and a grey cover for the lever bar.



Carospeed Indicator is only available with a black handle.



Product information

CAROSPEED	CAROSPEED E	CAROSPEED T	CAROSPEED TE
<ul style="list-style-type: none"> ● Black knob ● Grey cover ● Hill-holder ● Set function for cruise control (option) 	<ul style="list-style-type: none"> ● Black knob ● Grey cover ● Hill-holder ● Set function for cruise control (option) ● El. regulator 	<ul style="list-style-type: none"> ● 2 different handle designs ● Grey cover ● Hill-holder ● Set function for cruise control (option) 	<ul style="list-style-type: none"> ● 2 different handle designs ● Grey cover ● Hill-holder ● Set function for cruise control (option) ● El. regulator
<ul style="list-style-type: none"> ● Wood/black knob ● Grey cover ● Hill-holder ● Set function for cruise control (option) 	<ul style="list-style-type: none"> ● Wood/black knob ● Grey cover ● Hill-holder ● Set function for cruise control (option) ● El. regulator 	CAROSPEED INDICATOR	
		<ul style="list-style-type: none"> ● Black knob ● Grey cover ● Hill-holder 	<ul style="list-style-type: none"> ● Direction Indicator switch ● Set function for cruise control (option)
Fulfils following demands: ISO/TC22 N 2072, 74/60/EEC, 78/632/EEC, 95/54/EC, VVFS1994:5			

Accessories

STEERING DEVICE
Makes steering easier when using only one hand.

CRUISE CONTROL
If the original vehicle has a cruise control this can be connected to the Carospeed set function for cruise.

ARMREST
An armrest prevents the arm from being strained when driving with Carospeed.

UNIVERSAL BRAKE PEDAL COVER
For increased safety a removable brake pedal cover is available. It prevents feet from being trapped under the pedal during drive.

MOUNTING BRACKETS
The Carospeed mounting bracket is fixed in the front seat bolt holes. Vehicle Specific or Universal.

Please
find more
information on
www.autoadapt.com



Autoadapt AB, Hedeforsvägen 6,
S-443 61 Stenkullen, Sweden
Tel: +46 (0)302 558 20 Fax: +46 (0)302 558 29
E mail: contact@autoadapt.se www.autoadapt.com

Autoadapt UK Ltd, Unit 4, Windsor Industrial Estate,
Rupert Street, Aston, B7 4PR Birmingham, UK
Phone: 0121 3335170 Fax: 0121 3596194
E mail: contact@autoadapt.co.uk www.autoadapt.co.uk

