

Mini-Hoists 40kg & 80kg



Our range of hoists come with a wide selection of options that power in/out and various lifting attachments.

Quick, effortless loading of your scooter or wheelchair.

Key features

- Lifts mobility vehicles weighing up to 80kg
- Tailor made fitting kits to suit most vehicles
- Simple lifting brackets to suit most scooters and wheelchairs
- Saves VAT on vehicles when fitted from new (ask for details)
- The same sequence applies when loading electric wheelchairs



1

With the scooter positioned as shown, attach the bracket to the scooter.



2

Simply press the up button and the scooter will be lifted up clear of the vehicles bumper.



3

Once the scooter has reached the correct height, press the 'in' button and the scooter can be manoeuvred into the vehicle.



4

Press the down button to lower the scooter into position. Apply restraint if fitted.

Quick, effortless loading of your scooter or wheelchair.

Key Benefits

- Lifts wheelchairs and scooters complete, dependent on the vehicle
- Powered directly from your vehicle's own battery
- Overload protection is standard on all our mini-hoists
- Can be removed easily (including power in/out version) to leave rear of vehicle clear. Does not apply to mounting frame.



Hoist Arm Folding Mechanism



Fitting Kit



Restraint



"I have to write to congratulate everyone at Autochair that was involved in providing me with a wonderful stress free experience in the recent purchase and fitting of the Hoist in my Zafira. From the original call to you, to the finished item, most professional, I cannot thank you enough. The name "Autochair" is the first name I will remember and recommend when anyone asks about hoists"

James Pye,
Stoke on Trent

- Lifting brackets are simple to attach, even for people with limited dexterity or strength
- Can be fitted to most vehicles
- Transfers easily from vehicle to vehicle using tailored fitting kits
- Guaranteed for twelve months (extra guarantee and service agreement available)



Lifting Bracket



Motability accredited.
NEMEDA UK member.
European Mobility Group member.

Autochair have over 30 years experience in the design and manufacture of vehicle adaptations to assist people with disabilities and reduced mobility - **call us today to find your new freedom with Autochair: 0800 214 045**

- Wheelchair Hoists
- Scooter Hoists
- Swivel Seats
- People Lifts
- Wheelchair Carriers
- SMARTDRIVING CONTROLS

Autochair reserve the right to change product specifications or appearance without prior notice through product development and improvement. We aim to ensure our customers are aware of this at the time of ordering.