



Autoadapt develops and produces products to make your vehicle better suit your needs

January 2006

Autoadapt gives Challenged Kids new possibilities



TURNY
Orbit



CARONY



TURNOUT



CAROSPEED

Autoadapt's strength today is that we solve the problem of getting in and out of vehicles. Our products give people with disabilities the opportunity to be more independent, with a minimum of assistance when wanting to go somewhere by car.



The Swedish company Autoadapt is one of the world's leading companies within product development, manufacture and distribution of product aids for vehicle adaptation. For 20 years our mission has been to offer safe, functional and well designed solutions for the disabled.

Since 2003 the USA based Bruno Independent Living Aids has been a partner in Autoadapt AB. This means we can continue to focus on product development, expand on an international market and become an even stronger partner to work with. Today we have some 90 employees and with that are one of Europe's largest companies within our industry. Our products are currently exported to more than 30 countries.

Please find more information on www.autoadapt.com



Autoadapt AB, Hedeforsvägen 6,
S-443 61 Stenkullen, Sweden
Tel: +46 (0)302 558 20 Fax: +46 (0)302 558 29
E mail: contact@autoadapt.se www.autoadapt.com

Autoadapt UK Ltd, Unit 4, Windsor Industrial Estate,
Rupert Street, Aston, B7 4 Birmingham, UK
Phone: 0121 3335170 Fax: 0121 3596194
E mail: contact@autoadapt.co.uk www.autoadapt.co.uk



INVENTIONS FOR QUALITY OF LIFE

CARONY Kids

The chair that grows with your child!

CARONY Kids does away with having to lift the child and simplifies transport to and from the car. Dock the back of the wheelchair towards the vehicle swivel base and release the seat from the wheel unit. Slide the seat over from the wheelchair rails to the rails of the vehicle swivel base. When the seat has been transferred you just swivel the seat facing towards the front, fasten the seat belt and you are ready to go.



The CARONY Kids is combined with either the swivel base TURNOUT or the car seat lift TURNY on which rails for sliding the seat are mounted.

The seat cushion of CARONY Kids can be length adjusted to follow the growth of the child. The cushion and backrest can be washed and changed, if you wish. When the child has grown out of the chair, it can be converted to adult size with the help of a kit from Autoadapt. Accessories, as well as alternative wheel unit colours, are available as options.

CARONY Kids is available in following alternatives:

- with **12" rear wheels**, which is a transport wheelchair.
- with **24" driving wheels** for children who can drive the chair by themselves.
- with **CARONY Fixed** wheel unit which is of a standard, fixed height. Used in combination with the TURNY.



Carony Kids with the Turnout swivel



The seat depth is 40 cm and can be extended by 5 cm. Neck cushion and body support can be adjusted to give best possible support to the child as it grows.

Carony Kids combined with **TURNY**

TURNY -the seat lift for vehicles with higher floor level (MPVs and minibuses). It swivels out the seat and lowers it outside the vehicle.



CARE BELT system

- Harness system which provides comfort and several combination options!

Careva STABILIZING BELT

- can be fitted in all vehicles with seats where there is an opening between the seat and the back.
- can be fitted in both front and back seats (but not in fixed back seats where special adaptation is required).
- can be used in child safety car seats (forward-facing) that have an opening between the seat and the back.
- can be used in conjunction with a belt cushion.
- can be fitted in all wheelchairs with a fixed back and an opening between the seat and the backrest.



The Careva belt is CE- marked (which indicates product compliance with the EU medical aids directive).




The mobility-impaired often require extra support or restraint to ensure a positive sitting position in cars, wheelchairs, child car safety seats, etc. The Careva belt system comprises stabilizing belts that can be combined in different ways, thus providing an excellent complement to conventional belts in restraining the torso against the backrest (or in special vehicles that lack a three-point belt). This also helps to prevent injuries during violent braking.

NB! The belt system should always be used together with the vehicle standard seatbelt!



The Spinal Belt is the base on which all the stabilizing belts are fixed.



Shoulder Harness Thoracic Belt Pelvic T-belt Combi Harness