


GETTING SEATED CAROLOCK



GETTING SEATED

CAROLOCK






The **Carolock** is a user-friendly wheelchair securement system for your vehicle, guaranteeing safety and worry-free stowing of your mobility device. It ensures that your device remains securely in place throughout your journey. The **Carolock** is built to support both independent drivers and passengers with reduced mobility.

Adjustments can be made to fit the needs of different wheelchair types, whether you have a newer or older model, a manual, or a powered wheelchair.






Once you reach your destination, return to your wheelchair, and release the device from the **Carolock** using a release handle attached to its mechanical system. By pressing the release handle, your device will be released safely and smoothly.






FEATURES

-  The height can be tailored to suit individual requirements and different wheelchairs.
-  The users mobility device will always be securely positioned.
-  The mechanical securment system makes it effortless to position and lock your wheelchair.

SPECIFICATIONS

-  Designed for mobility devices up to 85kg.
-  Height - 122mm
-  Width - 88mm
-  Cycle Test Approved 173.000
-  Pull Test Approved

BENEFITS

-  Can be combined with the Turny 6 Way.
-  Flexibility to decide how and where the Carolock can be installed.
-  Get in your preferred seating position and concentrate on your travels without any safety concerns!

