


GETTING SEATED

TILDA



GETTING SEATED

TILDA

*The Tilda can be paired with the **Turny Evo** or **Turnout** swivel seats for maximum effectiveness.*

Tilda is an ergonomic tilt function with many uses. Installed under the car seat it allows the seat to be tilted in three steps, giving the user options in seating position for pressure relief and comfort when using a swivel seat. It can also be used to lower the user's head enough to clear the roof when swiveling in or out of the vehicle.





Tilda combines excellently with the seat swivel TURNOUT and makes it even easier to get in and out of the car. By keeping the feet on the footrest they automatically follow over the threshold.

This product improves the head space when sliding the seat between the Carony wheel unit and the Turnout or Turny swivel seat inside the vehicle.




The optional footrest can also be used to provide help to get the feet easily over the threshold when swiveling inside the vehicle. Tilda is also an excellent option for a safe and comfortable seating position.



FEATURES

-  Comes with an optional foot rest for comfort.
-  Is capable of combining with the Carony wheelchair base.
-  Three adjustable positions to suit your needs.
-  Can be used with the Turnout or the Turny Evo to increase knee space.

SPECIFICATIONS

-  Weight: 7 kg
-  Height: 2.4 cm
-  Tilting Angles: 0°, 7°, 11°, 15°