


GETTING SEATED TURNY MANUAL



GETTING SEATED

TURNY MANUAL







The **Turny Manual** is an easy to install system that is already fastened to a rail system, meaning the mechanical side of installation is extremely simple.




Part of the Turny family, this seat allows a user to conveniently transfer from outside to inside a vehicle with little to no hassle. With the use of a basic remote, the **Turny Manual** extends completely outside the vehicle, allowing occupants to move from wheelchair to seat, or even just into the vehicle without any struggle to access.

This system can be removed and installed into new vehicles as needed and can be modified to fit the driver and passenger side doors. Users are able to maintain their independence and eliminate any issues they have with getting in and out of their vehicle.



Features

-  Easy to Install
-  Slim Design
-  Comfort Adjustment
-  No Rails

Specifications

-  44mm more narrow than the Turnout model.
-  Rails are already built in.
-  Three positions for lengthwise adjustment.

Benefits

-  Simple to operate
-  Gives the user mobility independence

