

WHEELCHAIR HOIST—CAROLIFT 90



WHEELCHAIR HOIST—CAROLIFT 90



OVERVIEW


- The Carolift 90 is a stylish wheelchair hoist designed for lifting mobility devices up to 90 kg into the back of a vehicle. The lift can be folded down to allow use in vehicles with sloping rear doors.

FEATURES

- Lifts up to 90kg
- Powered lifting and swivel
- Remote control operation
- EMC and crash test approved

SPECIFICATIONS

- Lift Capacity:** 90kg
- Weight :** 15 (Including mounting bracket)
- Voltage:** 12V
- Current Consumption** : < 10A (at max load)



A = 106 cm including mounting base. Height adjusted by cutting

B = 109 cm. Length adjusted by cutting

C = Max 141 cm folded

D = Max 58 cm