


GETTING SEATED CARONY GO



GETTING SEATED CARONY GO



Transfer wheelchair system



The **Carony** is a system that transforms the car seat into a wheelchair and allows for a vehicle transfer that eliminates any manual handling of the occupant. The **Carony** base docks into the swivel seat in the vehicle, either the Turny or Turnout variants, and the entire seat easily slides from the car onto the wheelchair base.




The **Carony Go** is the full electric wheelchair version. This unit can also be used with an attendance drive for carer assisted transport when required.

There are 3 seating options for the wheelchair to suit the occupant's needs: **Bev seat, Compact seat or Recaro seat.**





This system is the alternative to modifying the vehicle for wheelchair access where the wheelchair occupant does not need a specialised wheelchair.



FEATURES



-  Left or Right applications
-  Wheelchair base can be adjusted in height electrically
-  Multiple options to cater to occupants seating and postural support

SPECIFICATIONS

-  Weight: 75 kg
-  Max user weight: 120 kg
-  Max Speed: 6.5km/h, Max Rang: 25kms
-  Crash tested and approved to ISO 7176 and ISO 10542 and CE marked



BENEFITS

-  Functionality of a regular wheelchair with comfort of an automotive seat
-  Suitable for both children and adults providing freedom of independence



 BraunAbility