

WHEELCHAIR HOIST—CAROLIFT 40



WHEELCHAIR HOIST—CAROLIFT 40

OVERVIEW

- The Carolift 40 provide excellent reach at the same time it is easy to operate through narrow opening due to it's flexible arm. It can fits into most vehicles due to it's small foot print.

FEATURES

- Lift up to 40kg
- Powered lifting and manual swivel
- Fits almost all car models
- Easy installation thanks to smart mounting base
- Operate via control on the hoist
- CE-marked

SPECIFICATIONS

- Lift Capacity:** 40kg
- Weight :** 11 kg (Including mounting bracket)
- Voltage:** 12V
- Current Consumption** : < 15A (at max load)



A = From cc pillar to lift belt, maximum — 795mm

B = From cc pillar to lift belt, minimum — 530mm

C = From hoist arm bottom edge to vehicle floor — 755mm

D = Total height , maximum — 875mm

Total height , minimum — 655mm

Length of lift belt — 1380mm

