


GETTING SEATED TURNY ORBIT



TURNY ORBIT



Programmable Lowering Swivel Seat System for transferring into your vehicle.



The **Turny Orbit** is a fully automatic, swivel seat where all movements are operated by the hand held remote.

The **Turny Orbit** allows the car seat to swivel out of the door opening and lower down to a height that is suitable for transfer from a wheelchair.

It is recommended to be used with one of our low-profile seats, so that you can optimise the head and knee room when inside the vehicle.

The Turny can also be used in combination with the Carony transfer wheelchair system.

FEATURES

- Fully powered swivel, lowering outside of the vehicle
- Suitable for all vehicles type such as SUV's, 4WD's, Vans and Campers
- Total user weight capacity of 150kgs

BENEFITS

- Reduces the effort required for vehicle transfers
- Provides you with the independence you deserve when travelling to and from your destination

SPECIFICATIONS

- Up to 480mm of vertical travel
- Easy to use back-up and emergency options
- Can be combines with Carony to eliminate any lifting during the transfer
- Optional footrest available for easier transfer



4/45 Salisbury Road
ASQUITH NSW 2077
(02) 9482 4572



www.mobilityengineering.com.au

sales@mobilityengineering.com.au