



GETTING SEATED TURNOUT



TURNOUT



Seat base that allows you to swivel out of the door for easy access






The **Turnout** is a swivel base that enables you to swivel the seat of the vehicle so that you are facing outside, making it easier to transfer in and out.

It is recommended to be used with one of our low-profile seats, so that you can optimise the head and knee room when inside the vehicle.







FEATURES

-  Manual and Electric versions available
-  2 or 4 door models to optimise swivel space
-  Total user weight capacity of 150kgs



SPECIFICATIONS

-  Easy to use back-up and emergency operations
-  Full backrest adjustability
-  Can be combined with Carony to eliminate any lifting during the transfer
-  Optional footrest available for easier transfer





Options for the **Turnout** include the Caroslide and Tilda, which can be used to further increase the swivelling room available in the vehicle.

Both *manual* and *electric* version is available.

The **Turnout** can also be used in combination with the Carony transfer wheelchair system.

BENEFITS

-  Reduces the effort required for vehicle transfers
-  Provides you with the independence you deserve when travelling to and from your destination



4/45 Salisbury Road
ASQUITH NSW 2077
(02) 9482 4572



www.mobilityengineering.com.au

sales@mobilityengineering.com.au